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Dear Customer

Thank you for purchasing a Dentaureum product. You have chosen a product that will provide you with consistent high quality, reliability and service.

Optimum performance of our dental technology products is based on working precisely according to the instructions issued by our Research and Development Department. We have therefore included a clear, detailed summary of each working stage in these instructions for use.

If you have any further questions, our team of orthodontic technicians will be pleased to offer advice.

Simply call:

Technical Advice,

fixed appliance technique Hotline Tel. no. +49 (0)7231/803-550

removable appliance technique Hotline Tel. no. +49 (0)7231/803-555

Information, instructions for use and safety data sheets for Dentaureum products are also available on the Internet at www.dentaureum.com.

Warranty

Dentaureum guarantees a faultless Quality of the products manufactured by us. These non-obligatory recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence in the processing on site. Eventually occurring claims for compensation may only be made up to the value of our goods.

Contents

	page
Introduction	14
Parts	14
Precise occlusal positioning of plaster models for horizontally divided appliances in the FKO-Split Post Fixator with the aid of the additional kit	15
Removal and trimming of models	16
Tips	16
Important notes	17
Setting the horizontal split level	17
Finishing the horizontally divided orthodontic appliance using the wire-spring activator after A.M. Schwarz as an example	18
Further working tips	19

The additional kit for the FKO-Split Post Fixator

Introduction

The additional kit for the FKO-Split Post Fixator is designed to aid in the efficient production of horizontally divided functional orthodontic appliances such as: the wire-spring activator after Schwarz, Siemons-Activator, Progenie-Activator after Wunderer, U-wire Activator after Karwetzky, Hansa-appliance after Hasund and many more. Due to the correct positioning of the lower model in the fixator, parallel to the occlusal surface, the additional kit aids in the simple production of anterior and lateral bite plates in the upper and lower jaw.

Parts

- Adjustment aid for the correct positioning of the lower model
- Metal mirror with fixation screw for producing full gloss, split horizontal planes
- Long cylinder with fixation screw to secure the position of the upper model in the Fixator (bite lock)
- Short cylinder with fixation screw to secure the position of the metal mirror in a bite lock in the intermaxillary area.
- Metal disc to compensate for the metal mirror thickness
- Complete instructions for use

Precise occlusal positioning of plaster models for horizontally divided appliances in the FKO-Split Post Fixator with the aid of the additional kit

1. The central screws on the fixator arms must be screwed in fully.
2. Remove the FKO-Split Post Fixator bite lock and replace with the long cylinder with fixation screws. Place the metal cylinder over the outer telescope, push down to the lower fixator arm and fix with the appropriate screws. (figs. 1 & 2).
3. Place the short cylinder with fixation screws onto the inner telescope, push down to the lower fixation arm and fix with the appropriate screws. (fig. 3).
4. Then place the metal mirror onto the inner telescope and fix with the appropriate screws. (fig 4).
5. Place the FKO-Split Post Fixator including the metal mirror into the adjustment aid (fig. 5).
6. Place the lower model onto the metal mirror so that the incisal edges of the front teeth and the disto-buccal cusps of the most posterior molars touch the metal surface (occlusal surface). If necessary fix the model with wax. (fig. 6).
7. Put the lower model into the FKO-Split Post Fixator and fix with articulating plaster. Ensure that the plaster on the lower model is flush with the lower fixation arm remove any excess plaster from the screw slits (fig. 7).
8. Once the plaster has fully set, take the adjustment aid out from the FKO-Split Post Fixator and remove the metal mirror (fig. 8).
9. Fix the upper model with the bite registration to the lower model (fig. 9).
10. Articulate the upper model in the FKO-Split Fixator in the usual manner (Fig.10).

Removal and trimming of models

1. Separate the upper fixator arm from the lower. Unscrew and remove the central screws (fig. 11).
2. Remove the models with their plaster bases from the fixator arms, taking care not to tilt to one side or the other (fig. 12).
3. Trim the plaster base right up to the plaster model (fig. 13).
4. When replacing the models, check that the fixator arms are clean (free of plaster and wax residues). If necessary, clean the fixator arms and model base using a steam cleaner.
5. Replace the models onto the fixator arms carefully without tilting.
6. Using a screwdriver, turn the central screws all the way back inside and so subsequently securing the models.
7. Replace the bite registration. Place the upper fixation arm onto the lower. Place the long metal cylinder onto the upper fixator arm and screw in tightly (fig. 14).

Tips

- seal the slits in the central screws to avoid plaster from entering.
- Spray the FKO-Split Post Fixator with a thin layer of silicon spray to allow easy removal of the models later.
- Any excess plaster on the sides of the upper and lower model can be easily removed using the plaster trimmer.
- Placing the plaster models and their bases into the archives saves time later. Repeating articulation is unnecessary should orthodontic appliance repairs be required. Simply slide the old models back onto the FKO-Split Post Fixator arms and secure with the central screws.

Important notes

- never strike the Fixator with a hammer or similar instrument in order to remove the models and plaster bases. This will distort the Fixator and damage its function.
- Never use metal instruments on the high gloss metal mirror (eg. metal plaster knives, metal spatulas etc.) in order to prevent unnecessary scratches. Damaged mirror surfaces create poor acrylic surfaces in the anterior and lateral occlusal region of the orthodontic appliance.

Setting the horizontal split level

1. Remove the bite registration
2. Place the metal disc, which compensates for the metal mirror, onto the outer round telescope between the long cylinder and the fixator arm. This guarantees that the finished orthodontic appliance has the same vertical dimension as the bite registration.
3. Positioning and fixing the metal mirror within the vertical centre of the intermaxillary bite lock. The position of the metal mirror is secured by lifting and fixing the small metal cylinder on the inner round telescope. (fig. 15).

Finishing the horizontally divided orthodontic appliance using the wire-spring activator after A.M. Schwarz as an example

1. Fix the expansion screw in the usual manner.

Caution:

- The direction of screw pressure must be absolutely parallel with the horizontally divided level (metal mirror / occlusal surface), in order to avoid uneven tilting when the appliance is activated.
- The acrylic screw adapter must not touch the metal mirror, shorten if necessary (fig. 16).

2. Bend and fix the wire elements in the usual manner (Fig.17).

Caution:

- The occlusal wires must not touch the metal mirror in any way (Fig.18,19).

3. Soak and apply a separating medium to the wires in the usual manner.

4. Apply the acrylic using either the dusting or dough method, in this case with the aid of a wax collar.

- Apply the acrylic in the lower appliance ensuring an over dimension of the occlusal level. (fig.20,21).
- Put on the metal mirror and adjust it as far as it will go, this will condense the acrylic in the lower occlusal area.(Fig.22)
- Apply the acrylic in the upper appliance with an over dimension of the occlusal level. (fig.23).
- Press the fixator together as far as it will go to the bite lock position, this will condense the acrylic in the upper occlusal area. (Fig.24)
- Polymerize the acrylic according to the manufacturer's instructions, with Orthocryl® 20 – 25 Minutes at 35° – 40°C and 2,2 bar pressure in the Polyclav® (Fig.25).

- After polymerization the occlusal level is totally flat and has a high gloss. It requires absolutely no correction. (Fig.26).
Finish and polish both acrylic bases in the usual manner taking care not to touch the occlusal surfaces of either the upper or lower appliance (Fig.27).
The finished and polished acrylic appliances: Fig.28,29,30.
- Bending, fixing and polymerization of the wire spring (Fig.31).
- The finished orthodontic appliance: Fig.32.

Further working tips

- A wax collar will prevent the flowing acrylic from leaking out.

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